

auroflex

with Aspheric Optic

CE
0197

“Your trusted
brand Auroflex is
now available
with
***Aspheric
Optic***”

Unfolding the Beauty of Vision . . . Again !



auroflex



A R A V I N D E Y E C A R E S Y S T E M

Raw material is imported from a renowned US supplier, who has been supplying blanks worldwide since 1998

An Aberration free optic:

- Auroflex Aspheric optic does not introduce any Spherical Aberration to the pseudophakic eye.
- It improves contrast sensitivity in dim and low light conditions and improves night driving vision.
- It does not induce other aberrations when decentered and does not reduce depth of focus to the patient.



Vision with normal spherical IOL



Vision with an Aspheric IOL

Specification :

Brand Name	Auroflex
Diopter range	10.0 D to 30.0 D (15.0D to 25.0D 0.5D increment)
Incision required	2.8 mm
Refractive Index	1.460
Haptic design	Dual Haptic
Optic Diameter	6 mm
Square edge	Yes, 360° posterior square edge design
Angulation	0°
A constant	118.0; 117.8 (IOL Master)
ACD	5.0 mm

Specifications comply with ISO 11979 series of standards

Other features:

- Dual haptic for better centration and adapts well to bear ocular bag contraction.
- Greater stability of an IOL due to large contact area
- Optimal refractive index to reduce glare
- Fenestration provided in the haptic enables fibrosis development for better fixation
- 360° square edge to reduce PCO
- Supplied with sterile disposable delivery system to ensure safety to the patient.

Information published in this catalogue is subject to change without notification

No 1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India

Phone: 91452 3096100 Fax: 91 452 2446200

E-mail: info@aurolab.com Web: www.aurolab.com