


SofPort™ AO

Advanced Optics Aspheric Lens System by Bausch & Lomb

MULTI-PIECE FOLDABLE ASPHERIC INTRAOCULAR LENS



OPTIC:

- 360° anterior/posterior square-edge
- 6.0 mm asymmetrical biconvex (1.3:1)
- UV absorbing
- Ophthalmic grade silicone (OGS2)
- Refractive index: 1.43

HAPTICS:

- Blue extruded PMMA
- Modified C loop
- 5° angulation

LENS CONSTANTS:†


- A-Constant: 118.0
- ACD: 5.0 mm
- Surgeon Factor: 1.22 mm

DIOPTERS:**

- 0 - 4.0D in whole diopters
- 5.0 - 30.0D in 0.5 diopters

LENGTH:

- 13.0 mm



LI61AO

* A-Constant, ACD and Surgeon Factor are estimates only. It is recommended that each surgeon develop his or her own values.

** 0 - 9.5D will be available in February, 2005.

† Exaggerated view of spherical and aspherical shapes.

Bausch & Lomb
SofPort™

Advanced Optics Aspheric Lens System

Start a resolution revolution with the new
Bausch & Lomb SofPort™ Advanced Optics
Aspheric Lens System.



CATARACT SURGICAL PRODUCTS

As a global leader in eye care, Bausch & Lomb is dedicated to perfecting vision and enhancing lives through innovative technology and design. To this end, we provide ophthalmic surgeons with a full array of products, including support for safe, minimally invasive cataract surgery.



IOL SYSTEMS

A wide range of advanced IOLs and inserter technologies including the new SofPort™ Advanced Optics Aspheric Lens System



VISCOELASTICS

A choice of viscoelastic products that deliver the right viscosity and clarity to satisfy different procedural demands, including Amvisc® Plus, Amvisc® and Ocucoat®



OPHTHALMIC INSTRUMENTS

The Storz® brand of precision-engineered microsurgical instruments, diamond knives and disposable surgical products designed to serve the needs of every surgeon for every patient



MICROSURGICAL SYSTEMS

When your technique meets advanced technology. The Millennium™ System delivers control through technological advancements in pump design and Dual Linear control of both aspiration and phaco power

To order Bausch & Lomb products, call 800-338-2020 or fax 800-362-7006.

180 East Via Verde Drive
San Dimas, CA 91773
www.bausch.com

®/™ denote trademarks of Bausch & Lomb Incorporated.
©2004 Bausch & Lomb Incorporated. 9/04 SU-3739

Bausch & Lomb
Perfecting Vision. Enhancing Life.™