



**SIDAPHARM<sup>®</sup>**  
Ophthalmic & Medical Disposable Products

---

# Ophthalmic Disposables

## SSpears

### PVA Eye Spears

- A soft, absorbing tool to drain excessive fluids
- Provided in a sterile pack of 5 pieces
- Ultrasmooth PVA material
- 5 year shelf life




#### Product Features:

- Lint free
- High density formulation
- Rapid fluid absorption
- Gamma radiation sterilized
- Double pouch packing

## SHooks

### Iris Retractors



Flexible hooks used to increase pupillary diameter during surgery. The attached silicon stopper allows better grip and easy adjustment of pupil size.

Provided in sterile pack of 5 pieces.



## SFluoro

### Fluorescein Sodium

Ophthalmic Strips  
1.0 mg



Uniform distribution ensuring optimal stain

Enhanced detection of:

- Corneal abrasions
- Foreign bodies
- Herpetic corneal ulcers

Individually packed, sterile ophthalmic strips



#### SIDAPHARM

21, Stageiriti & 24, Em. Fill Str.  
Thessaloniki, GR- 54352, Greece, European Union  
T: +30 2310 906660 info@sidapharm.gr  
F: +30 2310 989846 www.sidapharm.gr